CYBERMED INC., ONDEMAND3D TECHNOLOGY INC.

# OnDemand3D Fusion Technology

---

## White Paper

December 2009

---

● USA ● Republic of Korea  
[www.ondemand3d.com](http://www.ondemand3d.com)

Copyright(c) 1999-2009, Cybermed. All right reserved.

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
## What is the Mutual Information(MI)

OnDemand3D fusion technology uses a basic concept from information theory, mutual information (MI), or relative entropy, as a matching criterion. Using MI in image superimposition is first suggested for registration of 3D CT, MR and PET brain images of single object (Maes F, Collignon A, Vandermeulen D, Marchal G, Suetens P. Multimodality image registration by maximization of mutual information. IEEE Trans Med Imaging 1997;16:187-98).

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With  $H(A)$  and  $H(B)$  being the entropy of  $A$  and  $B$ , respectively,  $H(A,B)$  their joint entropy, and  $H(A|B)$  and  $H(B|A)$  the conditional entropy of  $A$  given  $B$  and of  $B$  given  $A$ , respectively.

The entropy  $H(A)$  is known to be a measure of the amount of uncertainty about the variable  $A$ , while  $H(A|B)$  is the amount of uncertainty left in  $A$  when knowing  $B$ . From (2), MI,  $I(A,B)$  is the reduction in uncertainty of the variable  $A$  by the knowledge of another variable  $B$ , or equivalently, the amount of information that  $B$  contains about  $A$ .

Let us suppose  $A$  to be the information of the cranial base in pre-treatment CBCT data, and  $B$  to be that in post-treatment data. As correlation between  $A$  and  $B$  is getting higher,  $H(A|B)$  or  $H(B|A)$  becomes smaller. If cranial base areas from two data are geometrically aligned, then  $H(A|B)$  or  $H(B|A)$  is minimized and MI is maximized in (2) and (3). We can produce superimposed CBCT data by maximizing MI with translation and rotation of images.

This method was also used by Cevidanes et al to obtain geometrical information from one software package that was subsequently used in another to compare pre-segmented surface models (article in Dentomaxillofac Radiol). We expanded the scope of this method to volume imaging and slice images and also refined the algorithm and user interface. This result is a very accurate and rapid method to superimposed various images one program.

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Cybermed Inc. Whitepaper
OnDemand3D Technology Inc.
Cybermed Inc.